



Detergent Measurements for One Gallon



DETERGENT SOLUTION 	ALCONOX PFS® Powered for Foodservice (1%)	ALCOJET® Low-Foaming Powdered Detergent (1%)	LIQUINOX® Critical Cleaning Liquid Detergent (1%)	CITRANOX® Acid Cleaner and Detergent (2%)	CITRAJET® Low-Foaming Liquid Acid Cleaner (2%)	*TERGAZYME® Enzyme-Active Powdered Detergent (1%)	DETONOX® Ultimate Precision Cleaner (3%)
ROUNDED TBSP 	2 tbsp	2 tbsp	2 tbsp	4 tbsp	4 tbsp	2 tbsp	6 tbsp
OUNCES 	1 oz	1 oz	1 oz	2 oz	2 oz	1 oz	3 oz

*For effective cleaning, only use warm water when making a Tergazyme solution. Do not exceed 130° F temperatures as it will degrade the enzyme.